THUASNE

LombaStab Dorso

| | | |
|----|--|----|
| fr | Ceinture de soutien thoraco-lombaire - «serrage facile et précis»* | 7 |
| en | Thoracolumbar support belt - "Easy, precise tightening" | 8 |
| de | BWS-/LWS-Stabilisierungsothese - „einfaches und präzises Justieren" | 10 |
| nl | Thoraco-lumbale rugbrace - "Eenvoudige en nauwkeurige aanspanning" | 11 |
| it | Cintura di sostegno toraco-lombare - "serraggio facile e preciso" | 13 |
| es | Faja de sujeción toracolumbar - "fijación fácil y precisa" | 14 |
| pt | Cinta de suporte toracolombar - "Aperto fácil e preciso" | 16 |
| da | Lændestøttebælte med thorakal støtte - "nem og præcis stramning" | 17 |
| fi | Rinta- ja lannerangan tukivyö - "Helposti ja täsmällisesti kiristettävä" | 19 |
| sv | Ryggstöd för bröstryggen - "enkel och precis åtdragning" | 20 |
| el | Ζώνη θωρακο-οσφυϊκής υποστήριξης - «εύκολη και ακριβής σύσφιξη»* | 22 |
| cs | Pás pro podporu hrudní a bederní páteře - „Jednoduché a přesné stažení" | 23 |
| pl | Orteza piersiowo-lędźwiowo-krzyżowa - „Łatwe i precyzyjne mocowanie" | 25 |
| lv | Krūšu-jostasvietas atbalsta josta - „Viegla un precīza savilkšana" | 26 |
| lt | Palaikomasis krūtinės ir juosmens diržas - „Lengvai ir tiksliai suveržiamas" | 28 |
| et | Rinna-nimpepiirkonna tugivöö - „Lihtne ja täpne pingutus" | 29 |
| sl | Oporni prsno-križni pas - »enostavno in natančno zapenjanje«* | 30 |
| sk | Hrudno-driekový podporný pás - „jednoduché a presné dotiahnutie" | 32 |
| hu | Háti-ágyéki fűző - „Egyszerű és pontos rögzítés" | 33 |
| bg | Торакo-лумбален колан - „Лесно и точно затягане" | 35 |
| ro | Centură de susținere toraco-lombară - „Strângere ușoară și precisă" | 36 |
| ru | Грудопоясничный ортез с системой простой и быстрой фиксации* | 38 |
| hr | Torakalno-lumbalni potporni pojas - „jednostavno i precizno stezanje" | 40 |
| zh | 调节简易准确的胸腰部支撑带* | 41 |
| ar | حزام لدعم الصدر وأسفل الظهر - "شد سهل ودقيق" | 43 |

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en

THORACOLUMBAR SUPPORT BELT - "EASY, PRECISE TIGHTENING"*

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Elastic fabric ①

Back panel with Coolmax lining

Adjustable back frame ②

Straps with hooks ③

Sternum support ④

2 abdominal panels with 2 finger-loops ⑤

Handles ⑥

Laces

Textile components:

Polyamide - Polyester - Elastane - Viscose - Polyoxymethylene - Polyethylene - Elastodiene

Non-textile components:

Aluminium - Polypropylene - SEBS - Polyamide - Polyurethane - Polypropylene PP

Properties/Mode of action

The device relieves low back pain, supporting the lumbar spine: compartmental lumbar unloading is therefore achieved via increased intra-abdominal pressure.

This device uses back pads, implementing a pulley-like system to distribute pressure.

The device maintains spinal hyperextension thanks to the adjustable straps, comfortable metal back frame and removable sternum support.

Indications

Stable vertebral fracture (traumatic and osteoporotic).

Scheuermann's disease.

Static disorders (kyphosis, lordosis).

Non-specific low back pain (acute, subacute and chronic phases).

Contraindications

Do not use in the event of an unstable vertebral fracture.

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom

an increase in arterial blood pressure is not recommended.

Precautions

It is recommended that a healthcare professional adjust back frame to the shape of the patient's back before first use.

If necessary, the healthcare professional who prescribed or delivered the device can cut the frame pouch and straps to the desired length with scissors.

Verify the product's integrity before every use.

Wear the device above thin clothing.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

Do not wear the product in direct contact with skin (risk of pinching).

Do not wear the device close to machines (risk of becoming caught).

Do not wear the product in a medical imaging machine.

Do not wear the product when driving a vehicle.

Do not wear the product while sleeping.

Consult a health professional if discomfort or pain occurs.

Undesirable side-effect

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Back frame adjustment

Measure the height of the patient's spine such that the lower edge of the back frame is level with the coccyx and the upper edge is level with thoracic vertebra T1 ①.

Loosen the four screws slightly and adjust the frame height before tightening the screws. ②③

Without using tools or heat source, adjust the frame to the curve of the patient's back.

Insert the fitted frame into the back pouch.

Shorten the pouch to the defined frame length by cutting it with scissors. ④

Strap adjustment

Ask the patient to tighten the straps until they are flat against his/ her body, such that they straighten the back slightly and stabilise the spine.

Shorten the straps used to tighten the back section using scissors. ⑤⑥

Fit the hooks back onto each strap.

Fitting the device

Open the belt and lay it out completely flat so that you can see the handles. ⑦

Check that the two tightening handles are fully retracted.

If necessary, unhook them and separate the two back parts by pulling the ends of the belt, then fasten the handles again.

Put on the brace in the same manner as a backpack. ⑧

Position the belt around the waist, placing your hands in the finger loops ⑨ located on the abdominal panels. ⑩

The tightening back pads must be positioned outwards (i.e. not against the skin).

The lacing section must be centred over the spine; the pads should be positioned on either side of the spine.

The bottom of the belt must stop at the top of the crease between the buttocks.

The Thuasne logo indicates which way up the belt should be put on.

Fasten the belt : right side over left side, without pulling too tight or applying any particular pressure. ⑪

Adjust the tightness using the two handles ⑫: pull them forwards then fold them flat against the front of the belt, where they attach like a touch fastener. ⑬⑭

These two handles can be used to adjust the level of compression as required throughout the day ⑮.

Position the sternum support ⑯ and adjust the tightness; it should be positioned over the sternum, above the chest, without compressing the thorax ⑰

Adjust strap tension ⑱ and place the hooks on the self-fastening surface. ⑲⑳

Avoid placing the hooks against the tightening system thread as they can cause damage.

Washing instructions

Remove the back frame before washing. Close the self-fastening tabs before washing. Machine washable at 30 °C. Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of the product with household waste (no special disposal conditions).

* Study conducted in-house on a panel of 13 people, July 2019.

This medical device is a regulated health product with CE marking.

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Κοντύνετε τη θήκη ανάλογα με το ορισμένο μήκος της ενίσχυσης κόβοντας τη θήκη με ψαλίδι. ©
Thuasne SAS – SIREN/RCS Nanterre 542 091 186 – capital 1 950 000 euros